



THANK YOU for ordering the Freedom Leg!

**Please visit our website for detailed instruction video on correct
adjustment fitting of your Freedom Leg at:**

www.freedomleg.com/fitting

Unlimited Phone Support at 1-888-816-8127

Freedom Leg Fitting and Adjustment



The calf strap is secured with buckle



The foot should be floating above the brace.

Shin pad can be adjusted vertically to allow for comfortable placement and to ensure toe will not touch the ground while walking.

To adjust the height, loosen bolts about 2-turns to allow nuts to disengage from indexing slots.



When the straps are adjusted properly, the hips should be level. This is important to be able to comfortably walk.

For the top two straps, adjusting the lengths on the back side of the leg will allow the hip to raise or lower to level the hips.

Additional Fitting and Usage information and videos can be found at www.freedomleg.com or call 1-888-816-8127

FWD MOBILITY

Patents Pending