



## How it works: THE FREEDOM LEG

*Wrap-around thigh cuff stabilizes and supports femur*

*Adjusting angle of thigh cuff allows hips to remain level ensuring proper symmetry and gait*

*Dual Carbon Fiber beams provide knee stability*

*Upper and lower straps work together to stabilize and support lower leg and can accommodate a boot or cast*

*Properly adjusted, the foot will float above the floor with zero weight from the knee down*

*External stirrup is compatible with boot or cast*

## No Limitations!

### **WALK WITHOUT CRUTCHES**

With your hands free you can accomplish all your normal day to day tasks without assistance.

### **ADDITIONAL KNEE STABILIZATION**

When additional knee stability is required the Freedom Leg is the answer.

### **SAFER ON STAIRS**

With your hands available to hold onto a handrail you can safely navigate stairs.

### **COMPATIBLE WITH CAST OR BOOT**

Unique external stirrup wraps around and under cast or boot.

### **CONTROL COUNTERFORCE ON LEG**

Adjusting the length of the brace give complete control of the weight on the injured leg.

### **BEFORE AND AFTER SURGERY**

Perfect solution both before surgery and after surgery.

### **LIGHTWEIGHT AND DURABLE**

Made from light weight and durable advanced composite materials. Total weight is 3 lbs.

## Who Can Use the Freedom Leg?

The simple answer is ANYONE with ankle, foot or knee injuries requiring non-weight bearing. Walk with zero weight on the injured lower leg.

- Femur Fracture
- Femoral condylar fracture
- Fusion of foot or ankle
- Proximal tibial fractures
- Midshaft tibial fractures
- Fractures of the foot and ankle
- Micro fracture of the knee
- Micro-fracture surgery
- Osteotomies
- Patellar fractures
- Post-op neuromas
- Achilles ruptures
- Plantar wounds
- Bunionectomy



### Keep Living!

Go beyond the limitations of traditional crutches or scooters. With the ability to accommodate a cast or boot the Freedom Leg gives medical professionals more options in the treatment of their patients. The off-loading on the injured leg can be adjusted to allow from 0 to 100 percent of the user's weight to meet the medical requirements of the patient.

888-816-8127  
www.freedomleg.com

## 100% Off-Loading 100% Hands Free

